

## A-level Chemistry exemplar for required practical 8

### Measuring the EMF of an electrochemical cell

#### Student sheet


#### Requirements

You are provided with the following:

- pieces of copper and zinc foil (about 2 cm × 5 cm)
- propanone
- 2.0 mol dm<sup>-3</sup> NaCl solution
- 1.0 mol dm<sup>-3</sup> CuSO<sub>4</sub> solution
- 1.0 mol dm<sup>-3</sup> ZnSO<sub>4</sub> solution
- emery paper or fine grade sandpaper
- two 100 cm<sup>3</sup> beakers
- plastic or glass U-tube
- cotton wool soaked in sodium chloride solution
- voltmeter (digital or high impedance)
- two electrical leads with connectors for the voltmeter at one end and crocodile clips at the other end
- samples of metals.

#### Suggested method for setting up a standard cell


- Clean a piece of copper and a piece of zinc using emery paper or fine grade sandpaper.
- Degrease the metal using some cotton wool and propanone.
- Place the copper into a 100 cm<sup>3</sup> beaker with about 50 cm<sup>3</sup> of 1 mol dm<sup>-3</sup> CuSO<sub>4</sub> solution.
- Place the zinc into a 100 cm<sup>3</sup> beaker with about 50 cm<sup>3</sup> of 1 mol dm<sup>-3</sup> ZnSO<sub>4</sub> solution.
- Lightly plug one end of the plastic tube with cotton wool and fill the tube with the solution of 2 mol dm<sup>-3</sup> NaCl provided.
- Plug the free end of the tube with cotton wool which has been soaked in sodium chloride. Join the two beakers with the inverted U-tube so that the plugged ends are in the separate beakers.
- Connect the Cu(s)|Cu<sup>2+</sup>(aq) and Zn(s)|Zn<sup>2+</sup>(aq) half-cells by connecting the metals (using the crocodile clips and leads provided) provided to the voltmeter and read off the voltage.



---

### Suggested method for measuring comparative electrode potentials of different metals

- Clean a piece of copper using emery paper or fine grade sandpaper.
- Connect the positive terminal of the voltmeter to the copper using a crocodile clip and one of the leads.
- Cut a piece of filter paper to about the same area as the copper, moisten the filter paper with the sodium chloride solution and place on top of the copper.
- Connect the second lead to the voltmeter and use the crocodile clip on the other end of the lead to grip a piece of another metal.
- Hold the metal against the filter paper and note the voltage reading and sign.



- Repeat steps (d) and (e) with different metals and record your results in a table.
- Write the conventional representation for each of the cells that you have constructed
- Suggest how you could construct the cell with the largest EMF from the metals provided.